

# IMPLANTS

## For Total Hip Replacement

# IMPORTANT PATIENT LEAFLET



**READ CAREFULLY THIS LEAFLET BEFORE SURGERY. IT CONTAINS  
IMPORTANT INFORMATION FOR YOU.**

### WHAT IS IN THIS LEAFLET

Intended Use.....	2
Devices Expected Lifetime .....	2
Contraindications.....	2
Materials Compositions And Distributions For Hip Implants.....	3
Implants you have been implanted with: .....	4
Specific warnings & MRI Information .....	6
Possible side effects .....	7
Warnings and precautions .....	7



## INTENDED USE

**The hip implant system is intended to replace the natural hip joint (total or hemi hip arthroplasty) through an artificial system.**

Patients with a mature skeleton may be implanted with a total/partial hip prosthesis in cases of severe joint pathology (cartilage wear, inflammation like arthritis or after an injury such as a bone fracture after a fall) and where other devices or treatments have failed.

Please ask your surgeon for details. The decision to implant a total/partial hip prosthesis is left to the surgeon after evaluation of the risk/benefit balance and after discussion with you.

---

## DEVICES EXPECTED LIFETIME

The expected survivorships for the following hip implant systems are at least at 10 years:

- 95% for primary total hip replacement with a standard stem,
- 88% for primary total hip replacement with a lateralized stem,
- 84% for revision total hip replacement with a standard stem,
- 94% for hemiarthroplasty.

It can be reduced or extended depending on your activity level or events that could compromise implants integrity (falls, accidents...).

**Please attend regular check-ups according to your surgeon's advice.**

---

## CONTRAINDICATIONS

At the time of the surgery, these conditions are not recommended:

- Pregnancy and breastfeeding (only for components made of CoCr alloy);
- Allergy to the implants materials;
- Infections;
- Severe mental, muscular, neurological or vascular deficiencies affecting the limb in question;
- Destruction of bone or poor bone quality, which may affect the stability of the implant.

**These conditions should be discussed with your surgeon. Any pathology (even if not listed above) must be mentioned to your surgeon beforehand as well.**

---

# MATERIALS COMPOSITIONS AND DISTRIBUTIONS FOR HIP IMPLANTS

Mass-to-mass ratio %

<b>Stainless steel</b> C: ≤ 0.08% / Si: ≤0.75% / Mn: 2.00% - 4.25% / Ni: 9.00% - 11.00% / Cr: 19.50% - 22.00% / Mo: 2.00% - 3.00% / Nb: 0.25% - 0.80% / S ≤0.01% / P: ≤0.025% / Cu: ≤0.25% / N: 0.25% - 0.50% / Fe: Balance <b>Hydroxyapatite coating &amp; titanium undercoating</b> Ca <sub>5</sub> (PO <sub>4</sub> ) <sub>3</sub> OH: 100% & C: ≤ 0.1% / H: ≤ 0.3% / Fe: ≤ 0.6% / N: ≤ 5.0% / O: ≤10.0% / Ti: Balance	<b>SATURNE</b> dual mobility cup - <b>Cementless</b> <b>SATURNE II</b> dual mobility cup - <b>Cementless</b>
<b>Stainless steel</b> C: ≤ 0.08% / Si: ≤0.75% / Mn: 2.00% - 4.25% / Ni: 9.00% - 11.00% / Cr: 19.50% - 22.00% / Mo: 2.00% - 3.00% / Nb: 0.25% - 0.80% / S ≤0.01% / P: ≤0.025% / Cu: ≤0.25% / N: 0.25% - 0.50% / Fe: Balance	<b>SATURNE</b> Tripolar cup - <b>Cemented</b> <b>E2</b> stem - <b>Cemented</b> <b>ACOR</b> monobloc stem - <b>Cemented</b> <b>EVOK</b> stem - <b>Cemented</b> <b>INITIALE</b> and <b>INITIALE</b> revision stems - <b>Cemented</b> <b>GENERIC</b> and <b>GENERIC</b> revision stems - <b>Cemented</b> <b>SPHERIC</b> bipolar cup Stainless Steel femoral head
<b>BIOLOX® Delta Ceramic</b> Al <sub>2</sub> O <sub>3</sub> : 60.0 - 90.0% / ZrO <sub>2</sub> + HfO <sub>2</sub> : 10.0 - 30.0% / HfO <sub>2</sub> in ZrO <sub>2</sub> : ≤5.0% Intended additives: ≤10.0% Total amount of impurities ≤0.2%	Ceramic femoral head Ceramic revision femoral head <b>HORIZON II</b> liner
<b>Cobalt-Chromium alloy</b> Cr: 26.00% - 30.00% / Mo: 5.00% - 7.00% / Fe: ≤0.75% / Mn: ≤1.00% / C: ≤0.14% / Ni: ≤1.00% / N: ≤ 0.25% / Co: Balance	Cobalt-Chromium femoral head
<b>Highly cross-linked Polyethylene</b> PE: 100% Possible traces of Ti, Ca, Cl, Al.	<b>AUSTRAL</b> liner
<b>Titanium alloy</b> Al: 5.50% - 6.75% / Va: 3.50% - 4.50% / Fe: ≤0.30% / O: ≤0.20% / C: ≤0.08% / N: ≤0.05% / H: ≤0.015% / Ti: Balance <b>Hydroxyapatite coating &amp; titanium undercoating</b> Ca <sub>5</sub> (PO <sub>4</sub> ) <sub>3</sub> OH: 100% & C: ≤ 0.1% / H: ≤ 0.3% / Fe: ≤ 0.6% / N: ≤ 5.0% / O: ≤10.0% / Ti: Balance	<b>HORIZON II</b> cup with holes - <b>Cementless</b>
<b>Titanium alloy</b> Al: 5.50% - 6.75% / Va: 3.50% - 4.50% / Fe: ≤0.30% / O: ≤0.20% / C: ≤0.08% / N: ≤0.05% / H: ≤0.015% / Ti: Balance <b>Unalloyed titanium porous coating</b> O: ≤0.40% / Fe: ≤ 0.50% / C: ≤ 0.08% / H: ≤ 0.05% / N: ≤ 0.05% / Ti: Balance	<b>AUSTRAL</b> cup - <b>Cementless</b>
<b>Titanium alloy</b> Al: 5.50% - 6.75% / Va: 3.50% - 4.50% / Fe: ≤0.30% / O: ≤0.20% / C: ≤0.08% / N: ≤0.05% / H: ≤0.015% / Ti: Balance <b>Hydroxyapatite coating</b> Ca <sub>5</sub> (PO <sub>4</sub> ) <sub>3</sub> OH: 100%	<b>EVOK</b> stem - <b>Cementless</b> <b>F.A.I.R.</b> stem - <b>Cementless</b> <b>ACOR</b> monobloc stem - <b>Cementless</b> <b>INTEGRALE</b> stem - <b>Cementless</b> <b>INTEGRALE</b> revision stem - <b>Cementless</b> <b>OPTIMAL</b> stem - <b>Cementless</b>
<b>Titanium alloy</b> Al: 5.50% - 6.75% / Va: 3.50% - 4.50% / Fe: ≤0.30% / O: ≤0.20% / C: ≤0.08% / N: ≤0.05% / H: ≤0.015% / Ti: Balance	Sleeve ( <i>used with Ceramic revision femoral head</i> ) Modular neck Threaded pin – Ø4,7 mm
<b>Polyethylene</b> PE: 100% Traces possibles de Ti, Ca, Cl, Al.	<b>SATURNE</b> liner for dual mobility cup



## IMPLANTS YOU HAVE BEEN IMPLANTED WITH:

*The surgeon needs to tick the boxes corresponding to the devices you have been implanted with.*

Associations with other devices have to be validated by Amplitude.

### ACETABULAR COMPONENTS

#### HORIZON II CUPS & LINERS RANGE

- ☐ **HORIZON II** Cup (with holes) – Cementless
- ☐ **HORIZON II** liner to be inserted



#### SATURNE CUPS & LINERS RANGE

- ☐ **SATURNE** dual mobility cup – Cementless
- ☐ **SATURNE II** dual mobility cup – Cementless
- ☐ **SATURNE** dual mobility cup – Cemented
- ☐ **SATURNE** liner for dual mobility cup to be inserted



#### AUSTRAL CUPS & LINERS RANGE

- ☐ **AUSTRAL** cup – Cementless
- ☐ **AUSTRAL** liner to be inserted



#### SPHERIC RANGE


- ☐ **SPHERIC** bipolar cup



# FEMORAL COMPONENTS

## FEMORAL HEADS RANGE

- ☐ BIOLOX® *Delta* Femoral head
- ☐ Metal (Stainless steel) Femoral head
- ☐ Metal (Cobalt-Chromium) Femoral head
- ☐ BIOLOX® *Delta* Femoral head with sleeve



## FEMORAL STEMS RANGE

- ☐ **INTEGRALE** femoral stem - Cementless
- ☐ **GENERIC** femoral stem - Cemented
- ☐ **INITIALE** femoral stem - Cemented
- ☐ **INITIALE** femoral stem - dysplasia - Cemented
- ☐ **E<sup>2</sup>** femoral stem - Cemented
  
- ☐ **ACOR** anatomic monobloc stem **standard** - Cementless
- ☐ **ACOR** anatomic monobloc stem **lateralized** - Cementless
- ☐ **ACOR** anatomic monobloc stem **standard** - Cemented
- ☐ **EVOK** femoral stem **standard collarless** - Cementless or Cemented
- ☐ **EVOK** femoral stem **standard collared** - Cementless
- ☐ **EVOK** femoral stem **lateralized collarless** - Cementless
- ☐ **EVOK** femoral stem **lateralized collared** - Cementless
- ☐ **EVOK** femoral stem **high offset collarless** - Cementless
  
- ☐ **F.A.I.R.** femoral stem **standard collarless** - Cementless
- ☐ **F.A.I.R.** femoral stem **standard collared** - Cementless
- ☐ **F.A.I.R.** femoral stem **lateralized collarless** - Cementless
- ☐ **F.A.I.R.** femoral stem **collared** - Cementless
  
- ☐ **INTEGRALE** revision femoral stem - Cementless
- ☐ **GENERIC** revision femoral stem - Cemented
- ☐ **INITIALE** revision femoral stem - Cemented
  
- ☐ **OPTIMAL** revision modular stem - Cementless (*with modular neck*)
- ☐ **OPTIMAL** reconstruction modular stem - Cementless (*with modular neck*)
- ☐ Threaded pin – Ø4,7 mm





## SPECIFIC WARNINGS & MRI INFORMATION

**The materials used in your implant may trigger security gates/scanners.**

**Please inform personnel about your implanted device.**

**The following technical information is intended for medical staff, to ensure that MRI equipment is used in the best possible conditions.**



Non-clinical testing demonstrated that the implants are MR Conditional. A patient with one of these devices can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 3T for a 15min-procedure.
- Maximum spatial gradient field of 1.500Gauss/cm.
- Normal operation mode only;
- The image artefact caused by the implants extends approximately 118.5mm from these implants. MR images may be blurred next to the implants.



### **MAXIMUM MEASURED TEMPERATURE RISE** with 3T MRI, for a 15 min procedure

<b>EVOK stems – cementless</b>	9.2°C (calorimetry WB-SAR of 4.71 W/kg)
<b>F.A.I.R. stems – cementless</b>	8.1°C (calorimetry WB-SAR of 4.08 W/kg)
<b>All other implants</b>	5.7°C (calorimetry WB-SAR of 2.06 W/kg)

<b>Specific weight restrictions for the</b> <b>Modular necks in association with OPTIMAL modular stems</b>	
<b>High L/M Neck</b>	<b>Lateralized neck version</b> Maximum patient weight: 90 kg
<b>L/M + 10.5 Neck</b>	<b>Lateralized and Medialized neck version</b> Maximum patient weight: 90 kg

<b>Specific weight restrictions:</b>		
<b>INITIALE stem</b> Dysplasia	<b>Size 121D</b> Maximum patient weight: 72 kg	
<b>ACOR stem</b> Lateralized	<b>Size 1</b> Maximum patient weight: 60 kg	
<b>EVOK stem</b> Standard	<b>Size 7</b> Maximum patient weight: 68 kg	<b>Size 8</b> Maximum patient weight: 85 kg
<b>EVOK stem</b> Lateralized	<b>Size 8</b> Maximum patient weight: 64 kg	<b>Size 9</b> Maximum patient weight: 89 kg



## POSSIBLE SIDE EFFECTS

You must be informed about the side effects in relation to the device that might occur. Those recorded by Amplitude are rare (inferior to 1 out of 1,000 cases). The following devices recorded occasional side effects (between 1 out of 1,000 cases and 1 out of 100 cases): EVOK stem<sup>1</sup>, BIOLOX® *Delta* femoral heads<sup>2</sup>, INTEGRALE/GENERIC stems<sup>3</sup>, INTEGRALE revision stem<sup>4</sup>, GENERIC revision stem<sup>5</sup>, INITIALE revision stem<sup>6</sup>, OPTIMAL<sup>8</sup>, ACOR monobloc cementless lateralized<sup>9</sup>, AUSTRAL<sup>10</sup>, M30NW and CoCr femoral heads<sup>11</sup>. The effects in **red** are recorded as 'frequent' for concerned devices (between 1 out of 100 cases and 1 out of 10 cases).

### **Please consult your surgeon if any doubt.**

- Allergy to implanted materials,
- Bone degeneration,
- Bone fracture<sup>1, 3, 4, 5, 6, 8, 10</sup>,
- Conflict between the components of the prosthesis,
- Difficult recovery after major surgery,
- Friction/conflict between the prosthesis and the surrounding tissue,
- Infection<sup>1, 2, 4, 5, 6, 8, 10</sup>,
- Neurovascular disorders,
- Noise from the prosthesis,
- Pain,
- Premature wear on joint surfaces,
- Prosthesis breakage,
- Prosthesis dislocation<sup>1, 4, 5, 10</sup>,
- Prosthesis displacement<sup>8</sup>,
- Prosthesis loosening<sup>4, 5, 9, 10</sup>,
- Tissue reactions to prosthetic debris and wear particles,
- Uneven legs length.

Other events linked to the surgery or current knowledge can also occur:

- Bone ossification abnormalities.
- Bruises.
- Phlebitis, fat or pulmonary embolism.
- Poor bone healing.
- Serious/ fatal complication of orthopaedic surgery, involving bone cement (only for cemented stems).

**You must also be aware of the risks and possible side effects of surgery - please ask your surgeon for further information.**

After a joint replacement, you should pay attention during your daily life: any blow, fall or accident may **jeopardize your recovery or your implant lifetime**. Please ask your surgeon for further information.

---

## WARNINGS AND PRECAUTIONS

You must be informed about the factors that could compromise the success of the surgery and post operative results:

- Overweight; Weight gain after surgery (especially for stems with weight restrictions).
  - Case history detailing infections and/or falls.
  - Metabolic disorders that reduce your resistance or induce progressive bone deterioration.
  - Local bone tumours.
  - Severe bone deformities.
  - Severe osteoporosis.
  - Addictive behaviour.
  - Playing sports intensively.
  - Playing risky sports or engaging in risky activities.
-

### **KEEP THIS LEAFLET**

You may need to read it again. If you have any further questions, ask your surgeon. If you experience any side effect, including **side effects not listed** in this leaflet, talk to your doctor or pharmacist. See page 7.

This leaflet has been written assuming the person receiving the implants will read it.

### **REPORTING ADVERSE EVENTS**

Any incident in relation to the device shall be reported to your surgeon, healthcare facilities or family doctor.



AMPLITUDE

11 cours Jacques Offenbach – ZA Mozart 2 – 26000 VALENCE – France

<https://implantcard.amplitude-ortho.com>



*BIOLOX® is a registered trademark of the Group CeramTec.*